
INSTRUCTION MANUAL

ELECTROPHORETIC BLOTTING SYSTEM

EBX-700

T A B L E O F C O N T E N T S

Important User Information	3-4	
Section 1: General Information		
1.1 Introduction and Features	5	
1.2 Specifications	5	
1.3 Safety	6	
Section 2: Description of Parts		
2.1 Unpacking and Components	7	
Section 3: Instructions for Use		
3.1 Blotting Unit Preparation	8	
3.2 Preparing the Gel for Transfer	8	
3.3 Preparation of Electroblothing Components	8-9	
3.3 Standard Electro-Blot Transfer	10	
3.3 Removing the Cassettes	10	
Section 4: Running Conditions		
4.1 Recommended Buffers, Membranes, & Running Conditions	11	
4.2 References	11	
Section 5: Maintenance of Equipment		
5.1 Care and Handling	12	
5.2 Maintenance	12	
Section 6: Ordering Information		
6.1 Ordering Information for EBX-700 Electrophoretic Blotting System	13	
6.2 Ordering Information for Related Products	14	
Section 7: Contact Information		15

IMPORTANT USER INFORMATION

This Instruction Manual will explain how to use this product safely and effectively. Please read and carefully follow the instruction manual in its entirety.

The triangle/exclamation mark symbol alerts the user of the product to important operational, maintenance, and/or warranty requirements.

The triangle/lightning bolt symbol alerts the user of the product to potentially hazardous electrical exposure.

Failure to adhere to the instructions could result in personal and/or laboratory hazards, as well as invalidate any warranty. Always turn off the DC power source prior to disconnecting power cords from the product. Disconnect power cords from the power source first and then from the product. For maximum safety, always operate this system in an isolated, low traffic area, not accessible to unauthorized personnel. Never operate damaged or leaking equipment.

WARRANTY AND LIABILITY

This product was produced utilizing the highest practical standards of materials, workmanship, and design. C.B.S. SCIENTIFIC warrants that the product has been tested and will meet or exceed published specifications. This warranty is valid only if the product has been operated and maintained according to the instructions provided.

C.B.S. SCIENTIFIC warrants this product to be free from defects in materials and workmanship under normal service for one year from date of shipment. If the product proves defective during this period, C.B.S. SCIENTIFIC will repair or replace it at our option, free of charge, if returned to us postage prepaid. This warranty does not cover: damage in transit, damage caused by carelessness, misuse or neglect, normal wear through frequent use, damage caused by solvent corrosion, damage caused by improper handling or user alteration, nor unsatisfactory performance as a result of conditions beyond our control. C.B.S. SCIENTIFIC shall in no event be liable for incidental nor consequential damages, including without limitation, lost profits, loss of income, loss of business opportunities, loss of use and other related damages, however caused, nor any damage arising from the incorrect use of the product.

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General Information

1.1 Introduction

C.B.S. Scientific now offers the EBX-700 4-Place Blotting System for performing electro-blotting. This unit provides the capability of blotting four gels simultaneously under identical temperature controlled buffer conditions. Units include 4 freezer blocks, and 4 blotting cassettes with 4 sponge pads.

Table 1: Features of EBX-700 Electrophoretic Blotting System

Model #	# of gels	Plate Dimensions (w x h)	Compatible Precast Gels	Cooling
EBX-700	4	Fits any gel up to 10cm x 10cm	Cambrex ClearPAGE Invitrogen	2 freezer blocks with stirring

1.2 Specifications

Constructions:

Buffer chamber, safety cover
electrode panels
Electrodes
Power cords
Safety Certification

Polycarbonate
 Polycarbonate
 Pure Platinum wire .012" diameter
 FR Urethane rated 7500VDC, 200mA, 65°C
 EN61010-1-1993 (IEC1010-1)

Table 2: Specifications

Model #	EBX-700
Shipping Weight	5 lbs
Overall Size	14(l) x 13(w) x 16(h) cm
Maximum Gel Size	10cm x 10cm
Recommended buffer volume	<ul style="list-style-type: none"> • 1.2 liters with no freezer blocks • 850 mls with 2 freezer blocks
Distance between electrodes	9 cm

1.3 Safety

Power to the EBX-700 Electrophoresis System is to be supplied by an external DC voltage power supply that must be ground isolated so that the DC voltage output floats with respect to ground. For any power supply used, the maximum specified operating parameters for the units are:

Maximum Limits

100 VDC
150 watts power
1500mA current
50°C ambient temperature

Current to the unit, provided from the external power supply, must enter the unit through the lid assembly, providing a safety interlock to the user. Current to the unit is broken when the lid is removed. **Do not attempt to use the unit without the safety lid, and always turn the power supply off before removing the lid, or when working with the unit in any way. Follow safety precautions specified by the power supply manufacturer.**

SECTION 2


Description of Parts

2.1 Unpacking and Components

Please verify that your unit comes complete with the following components:

EBX-700 4-Place Electrophoretic Blotting System:

- Blotting Chamber with 4 gel capacity
- Safety cover with attached DC power leads
- 4 freezer blocks
- 4 extra sponge pads (not shown)
- 4 blotting cassettes with sponge pads (not shown)

Instructions for Use

3.1 Blotting Unit Preparation


1. Place 4 freezer blocks in freezer and allow enough time to freeze.
2. Place the blotting chamber on a level work surface in an authorized work area.
3. Consult applications table 4.1 for buffer system and power settings.
4. Place, do NOT drop, a magnetic spin bar (not supplied by C.B.S.) in the bottom of the tank. This maintains pH, buffer and heat circulation during the electroblot.
5. Place tank on top of the magnetic stirrer.
6. Buffer should be pre-cooled to 10°C.

3.2 Preparing the Gel for Transfer


The EBX-700 can hold up to 4 gel cassettes simultaneously with space provided for a magnetic stir bar to allow buffer circulation and heat exchange for uniform transfers. When the resolving gel electrophoresis is complete, proceed with staining and photo-documentation if applicable. Cut one corner off the gel so that correct orientation is maintained throughout the procedure.

3.3 Preparation of Electroblotting Components

1. CBS membranes are supplied with the blotting paper cut to the correct size for ClearPAGE gel transfers. If not using CBS membranes, cut transfer membrane to size of gel or sufficient dimensions to cover the relevant bands. Also, pre-cut Whatman-type 3MM filter paper to same size and soak in transfer buffer until completely saturated (15 to 30 minutes).
2. Pre-wet membrane chosen according to manufacturer's instructions. Nylon or Nylon-supported nitrocellulose should be soaked in ddH₂O. PVDF in MeOH. Equilibrate all types in transfer buffer for 3 minutes. Submerge a single foam pad and gently tease to release trapped air bubbles until pad is completely saturated.
3. Pre-equilibrate gel in 1x transfer buffer to be used prior to electroblotting for 5 to 30 minutes depending on gel thickness.
4. Open the transfer cassette as shown below and lay it flat on the bench.

5. Assemble blotting stack as shown below. With cassette wide open assemble components on black side in the following order: buffer saturated sponge pad, gel equilibrated in transfer buffer, buffer saturated transfer membrane, then buffer saturated blotting paper. Smooth with gloved finger or roll with glass rod to be sure no bubbles exist between the gel and the transfer membrane.

6. **Insert blotting cassettes into the lower reservoir. Red side of cassette should face toward red electrode panel.**

7. Place 2 freezer blocks into the space provided in the lower reservoir. Depending on transfer conditions buffer temperature rise can be limited by exchanging freezer blocks half way through run.

3.4 Standard Electro-Blot Transfer

1. Fill the chamber with cold transfer buffer. (Buffers should be prepared fresh with reagent grade chemicals and pre-cooled). The transfer buffer should come up to the top of the platinum labyrinth. The buffer should not come in contact with the banana plugs when the gel cassette sandwich(es) are immersed in the unit.
2. Align safety cover over the unit and carefully attach.
3. Connect the leads to the power supply, matching the color-coded (+) red to red and (-) black to black. **See Section 4.1 for recommended power conditions.**
4. Transfer times will vary according to several parameters. Optimization of electro-blotting transfers must be determined empirically. Keep in mind the following principles that govern the movement of molecules of gel electrophoresis:
 - Thicker or higher percentage gels will take longer to transfer than thinner or lower percentage gels.
 - Large molecules will need extended transfer times to completely transfer.
 - Actual transfer times for defined conditions can be approximated by running molecular weight standards.

3.5 Removing the Cassettes

1. Turn the power supply off and disconnect the leads from the power supply.
2. Remove the safety cover from the unit, by placing thumbs on clear posts next to red & black connectors, then pushing down while pulling up with fingers under lid. **DO NOT pull on power cords.**
3. Gently lift the cassette from the unit. **Always wear gloves, eye protection and protective clothing** if buffer and/or gel contain Ethidium Bromide. Ethidium Bromide is a powerful mutagen; gloves, eye protection and protective clothing should always be worn when handling the gel or buffer solutions. See Material Data Safety Sheets.
4. Mark the orientation of the membrane with a pencil or by cutting off a corner and take apart the cassette carefully.
5. Process membrane according to type of transfer and manufacturer's recommendations.

SECTION 4

Applications & Running Conditions for EBX-700 Electro-Blotting

4.1 Recommended Buffers, Power Settings and Transfer Times

1. Using 2 freezer blocks and chilled buffer, start run at 60Volts(c) @ 370mA for 1-3 hours. Run time will vary depending on thickness of gel and sample size.
2. Buffers: ClearPAGE Transfer Buffer (cat. # CB82500) made up as follows:

COMPONENT	AMOUNT NEEDED FOR EBX-700
ClearPAGE Transfer Buffer	100 ml
Methanol	200 ml
Ultrapure water	720 ml

If you prefer to make up your own transfer buffer please make up 1X Towbin Buffer plus 0.1% SDS:

Protein Electro-Blotting Buffer:

TG-SDS (1X) 20% MeOH:

0.025M Tris base
0.192M Glycine
0.05-0.1% (w/v) SDS
20% (v/v) Methanol

4.2 References

- Ausubel, F.M., Brent, R., Kingston, R.E., Moore, D.D., Seidman, J.G., Smith, J.A., Struhl, K. (ed.) (1993). *Current Protocols in Molecular Biology*. Vol. 2, Greene Publishing Associates, Inc. and John Wiley & Sons, Inc., Ch.10.
- Burnette, W.N. (1981). Western blotting: Electrophoretic transfer of proteins from sodium dodecyl sulfate-polyacrylamide gels to unmodified nitrocellulose and radiographic detection with antibody and radioiodinated protein A. *Anal. Biochem.* 112:195-203.
- Peluso, R.W. and Rosenberg, G.H. (1987). Quantitative electrotransfer of proteins from sodium dodecyl sulfate polyacrylamide gels onto positively charged nylon membranes. *Anal. Biochem.* 162:389-398.
- Perides, G. Plagens, U., and Traub, P. (1986). Protein transfer from fixed, stained and dried polyacrylamide gels and immunoblot with protein A-gold. *Anal. Biochem.* 152:94-99.
- Tesfaigzi, J., Smith-Harrison, W., and Carlson, D.M. (1994). A simple method for reusing western blots on PVDF membranes. *Biotechniques.* 17:268-269.
- Towbin, H., Staehelin, T., and Gordon, J. (1979). Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: Procedure and some applications. *Proc. Natl. Acad. Sci. U.S.A.* 76: 4350-4354.
- Sambrook, J., Fritsch, E.F., Maniatis, T. (1989). *Molecular Cloning. A Laboratory Manual*. 2nd ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York. 18.47-18.61.

SECTION 5 Maintenance of Equipment

5.1 Care and Handling

The plastic components of the EBX-700 Electrophoretic Blotting System are fabricated from polycarbonate. Electrodes and connectors are made from pure platinum, stainless steel, and nickel plated brass. As with any laboratory instrument, adequate care ensures consistent and reliable performance.

After each use, rinse all parts with de-ionized water. Wipe dry with a soft cloth or paper towel, or allow to air dry. Whenever necessary, all components may be washed gently with water and a non-abrasive detergent, and rinsed and dried as above. *Never* use abrasive cleaners, glass cleaning sprays or scouring pads to clean the components, as these will damage the unit and components.

Additional precautions:

- Do not autoclave or dry-heat sterilize the apparatus or components.
- Do not expose the apparatus or components to phenol, acetone, benzene, halogenated hydrocarbon solvents or alcohols.
- Avoid prolonged exposure of the apparatus or components to UV light.
- Do NOT treat with diethylpyrocarbonate (DEPC)-treated water for extended periods at 37°C. A brief rinse with DEPC-water is sufficient after a thorough wash.

5.2 Maintenance

The following inspection and maintenance procedures will help maintain the safety and reliable performance of the EBX-700 Electrophoretic Blotting System. Replacement parts can be ordered by calling 1-800-243-4959 or by contacting your local distributor.

- Banana plugs and power cords should be inspected regularly. If the banana plugs become loose or do not feel friction tight replace the plugs or power cords.
- Should power cord assemblies (connectors, wire or shrouds) show any signs of wear or damage (e.g. cracks, nicks, abrasions, or melted insulation), replace them immediately.
- The platinum wire is secured to the banana jack by compression between a stainless washer and the jack nut. The nut/washer interface should be tight and free of corrosion.

6.1 Ordering Information

Part number	Description
EBX-700	4-Place Blotter, CE. Kit includes: lower reservoir, safety cover with attached leads, 4 blotting cassettes, 4 sponge pads, 4 freezer blocks, anode, cathode, and instruction manual
Accessories	
EBX-BC-700	Accessory Blotting Cassette, includes sponge pad
EBX-SP-700	Accessory Sponge Pads (set of 4)
EBX-FB-700	Accessory Freezer Block
HP19001	ClearPAGE Blot Sandwich NC 90x85mm 2/pk. 0.2 μ Nitrocellulose/1mm blotting paper
HP19020	ClearPAGE Blot Sandwich NC 90x85mm 20/pk. 0.2 μ Nitrocellulose/1mm blotting paper
HP29301	ClearPAGE Blot Sandwich PVDF 90x85mm 2/pk. 0.2 μ PVDF/1mm blotting paper
HP29320	ClearPAGE Blot Sandwich PVDF 90x85mm 20/pk. 0.2 μ PVDF/1mm blotting paper
EBX-EA-700	Replacement Anode
EBX-EC-700	Replacement Cathode
Power Supplies	
EPS-300-II	Mini Power Supply with timer, CV or CC, 10-300V, 110V/60Hz, current range: 4-500mA, 90 watts
EPS-300-IIV	Mini Power Supply with timer, CV or CC, 10-300V, 220V/50Hz, current range: 4-500mA, 90 watts

6.2 Ordering Information for Related Products

Part number	Description
DCX-700	Dual Cool System, CE. Fits precast gels or glass plate dimensions of 10x10 or 10x9cm(h). Kit includes: lower reservoir, safety cover with attached leads, core, 2 blotting cassettes with sponge pads, 2 freezer blocks, single gel adapter plate and instruction manual
Accessories	
EBX-BC-700	Accessory Blotting Cassette, includes sponge pad
EBX-SP-700	Accessory Sponge Pads (set of 4)
EBX-FB-700	Accessory Freezer Block
DCX-SP-109	Accessory Shim Plate for Cambrex 10x9 gel
DCX-AP-1010	Additional Single Gel Adapter Plate
Combs	
MVX-0701	Comb, 0.75mm x 1 well
MVX-1001	Comb, 1.0mm x 1 well
MVX-1501	Comb, 1.5mm x 1 well
MVX-0710	Comb, 0.75mm x 10 well
MVX-1010	Comb, 1.0mm x 10 well
MVX-1510	Comb, 1.5mm x 10 well
MVX-0714	Comb, 0.75mm x 14 well
MVX-1014	Comb, 1.0mm x 14 well
MVX-1514	Comb, 1.5mm x 14 well
Spacers	
DCX-S7510R	Gel Wrap Spacer Set, 0.75mm
DCX-S1010R	Gel Wrap Spacer Set, 1.0mm
DCX-S1510R	Gel Wrap Spacer Set, 1.5mm
Gel Wrap	
DCX-E7510	Gel Wrap Gasket, 0.75mm
DCX-E1010	Gel Wrap Gasket, 1.0mm
DCX-E1510	Gel Wrap Gasket, 1.5mm
Glass Plates	
DCX-P1010R	Gel Wrap Glass Plate Set, 10cm x 10cm

CALL FOR INFORMATION ABOUT OUR NEW LINE OF PRECAST GELS & ACCESSORIES....

ClearPAGE™ *Precast Gels
& Accessories*

NOTES

CONTACT US

Telephone:
Local or International
858-755-4959
Toll Free: 858-755-0733

Fax: 858-755-0733

Check out our new and improved online ordering!

www.cbsscientific.com

E-mail address:
sales@cbssci.com

Mailing address:
C.B.S. Scientific Company
P.O. Box 856
Del Mar, CA 92014

Shipping address:
C.B.S. Scientific Company
10805 Vista Sorrento Pkwy
Suite 100
Solana Beach, CA 92121

Credit Card Options:
Visa/Mastercard
Discover/American Express