

**PURIFYING TOTAL RNA FROM PLANT TISSUE OR FILAMENTOUS FUNGI****1 Sample homogenization**

See manual for recommended sample amounts and homogenization methods.

**2 Cell lysis**

Add 350µl Lysis Buffer RLY and 3.5µl β-ME to a maximum of 100mg ground tissue and vortex vigorously.

If lysate solidifies when adding Lysis Buffer RLY, use 350µl Lysis Buffer RLS instead.

**3 Filter lysate**

Place ISOLATE II Filter (violet) in a 2ml Collection Tube (supplied).

Load lysate and centrifuge 1 min at 11,000 x g.

Transfer filtrate to a new 1.5ml microcentrifuge tube (not supplied).

If visible pellet forms, transfer supernatant avoiding any pellet to a new 1.5ml microcentrifuge tube (not supplied).

**4 Adjust RNA binding conditions**

Discard ISOLATE II Filter and add 350µl ethanol (70%) to homogenized lysate.

Mix by pipetting up and down 5 times.

Alternatively, transfer flow-through into a new 1.5ml microcentrifuge tube (not supplied), add 350µl ethanol (70%) and mix by vortexing (2 x 5s).

**5 Bind RNA**

Place ISOLATE II RNA Plant Column (blue) in a 2ml Collection Tube.

Load lysate onto column and centrifuge 30s at 11,000 x g.

Place column in a new 2ml Collection Tube.

**6 Desalt silica membrane**

Add 350µl Membrane Desalting Buffer (MEM).

Centrifuge at 11,000 x g for 1 min to dry membrane.

**7 Digest DNA**

Add 10µl reconstituted DNase I to 90µl Reaction Buffer for DNase I (RDN).

Mix by gently flicking tube.

Apply 95µl DNase I reaction mixture directly onto center of silica membrane.

Incubate at room temperature for 15 min.

## 8 Wash and dry silica membrane

### 1<sup>st</sup> Wash

- Add 200µl Wash Buffer RW1.  
Centrifuge 30s at 11,000 x g.  
Place column into a new 2ml Collection Tube.

### 2<sup>nd</sup> Wash

- Add 600µl Wash Buffer RW2.  
Centrifuge 30s at 11,000 x g.  
Discard flow-through and place column back into Collection Tube.

### 3<sup>rd</sup> Wash

- Add 250µl Wash Buffer RW2.  
Centrifuge 2 min at 11,000 x g to dry membrane completely.  
Place column into a nuclease-free 1.5ml Collection Tube (supplied).

## 9 Elute RNA

Add 60µl RNase-free water (supplied) directly onto center of silica membrane.  
Centrifuge at 11,000 x g for 1 min.